Cell biophysics practical: Migration of cells in confined channels

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Day 1: Experiments

Tag 1: Messungen

Foreword: Along the description of the experiments you will find in bold some questions. As preparation we ask you to think about those questions and try to answer them. Your answers will also appear in your report. Along your answers we ask you to provide the source of your reflections (book title and authors, websites, lectures...).

During the experiments you should take notes of what you do since some things might differ from the protocol. This written protocol of yours will be looked at after the FoPra.

You should come some days before your appointment to watch an introductory video presenting the experiments. They are especially explaining those experiments for which you won't have the time to do them yourself and which we have done for you beforehand.

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The goal of our experiments is to observe, characterize and compare cells migrating in confined environment.

Cell migration is a way for a cell to move from a place to another. This is particularly important in many biological processes such as embryonic development, wound healing or immunity. Depending on their types, cells can move in different manners (migration mode). The two main modes studied for single cells are the mesenchymal migration and the amoeboid one.

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For mesenchymal migration, cells need to adhere to their environment (the working surface in the case of 2D migration, or the surrounding matrix in 3D). In short, the front of the cell adhere more firmly while the rear contracts. The force produced by the rear is then transmitted to the substrate via the adhesion points, producing a net movement of the cell.

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The case of amoeboid migration is more complex and its mechanisms are not yet fully understood. Nevertheless it is clear that for this migration to occur, cells don't need to adhere and the cytoskeleton needs to be contracted at the rear of the cell. One model for this migration is to see the cell as chimney climbers (Hawkins *et al.* PRL 102 (2009)). The rear contraction produces an increase of pressure within the cell that will be transmitted as a normal force to the surroundings. The pressure from the rear is then transformed as a flow towards the front, allowing the cell to move.

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One key component of both migration modes is the need to produce a force at the rear of the cells. In both cases it is achieved by a contractile force generated by the motor protein myosin II on the cytoskeleton.

The goal of this experiment will be to assess the importance of the myosin II in the case of amoeboid migration.

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To achieve this we will use a new setup composed of microfabricated channels of 5 microns size (Faure-André et al. Science 322 (2008)). The molds for the channels are built through microfabrication. This step is beyond the scope of this practical and you will be given already made ones for your experiments.

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Find or draw a simple scheme presenting the principles of Soft-photolithography.

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Biophysics studies often make use of PDMS chips. PDMS is an elastomer that is easy to prepare in laboratories. It is composed of a base polymer and a crosslinker which are mixed together. This mix will harden by heat.

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
From a material point of view, what is an elastomer. Specifically for PDMS find its advantages for biophysics that makes it so widely used in biophysics.

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For these experiments the microsystem used is presented on Scheme 1. This chip allows the visualization of cells migrating in short channels of 800 μm of length, and a square section of ca 5 by 5 μm . To ease the cells entry, there are two loading channels of 50 μm height on both ends of the migrating channels. Finally, side channels exit out of the chip to help cell medium to diffuse within the channels.

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Scheme 1 : Drawing of the PDMS chip

Schema 1: Zeichnung des PDMS-Chips

Your first task will be to create such a chip and clean it from the PDMS pieces you will be given at the beginning of the day. This corresponds to the first section of the protocol. At the end of that step you will have a PDMS piece with channels that are semi open. The next step will be then to close the channels binding the chip to a glass surface.

One advantage of the PDMS is that its surface can be easily modified. By using this property it is possible to bind it covalently to glass by plasma treatment.

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Recall the definition of a plasma. What will happen to PDMS and glass in plasma? How can they covalently bind? Best would be to provide a scheme of this process.

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In the second section of the protocol you will find the procedure to bind PDMS to glass with plasma treatment.

Before injecting the cells, the chip needs to be functionalized. Thanks to its known chemistry, PDMS can be functionalized with numerous different molecules. It is then possible for the cells to be in contact with adherent or repellent cues.

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For these experiments the channels will be coated with fibronectin, a protein that specifically interacts with different integrins to form focal adhesions and promote cell adhesion. While unnecessary for amoeboid cell migration inside the channels, this coating will help the cell to reach the channels' entries. The procedure is detailed in the section 3 of the protocol.

Once the chip is ready, cells can be loaded and their migration recorded. For these experiments, you will work with activated HL60 cells. HL60 is a human cell line derived from acute promyelocytic leukemia cells. Non-activated cells are proliferating in suspension and become adherent only after activation/differentiation. These differentiated cells (similar to neutrophils=white blood cells) are able to migrate.

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What is a cell? What is called a cell line? How do we culture cells (requirements, safety...)?

Cite at least 3 main organelles of a cell and their purpose.

What is called cytoskeleton and what are the main protein that compose it?

What are molecular motors, what's the role of Myosin II ?

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In order to record the migration of the cells in the channels, Höchst fluorescent staining will be used. This molecule binds to DNA which means it will specifically labels cell nuclei.

Once labeled the cells will be recovered from the cell culture dish and counted. They will be splitted in two sets of cells: one control condition and a second condition containing the drug Y27632 which interferes with the binding of myosin II to actin. The activity of the motor protein is hence reduced. Finally, cells for both conditions will be loaded into the channels. These steps correspond to the fourth section of the protocol.

The last step of the experiment is the recording of the migration by microscopy. To that purpose, we will use time-lapse, multi-position, and fluorescent microscopy. The procedure won't be detailed in the protocol part and will be explained during the experiments together with us.

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Recall the main optical components of a microscope. What is called inverted microscopy?

What is the principle of fluorescence? For fluorescent microscopy in cells, what are the main problems and what solution do we have to solve them? (and no, the answer to that one is not in the German Wikipedia ;-)) - but feel free to add it there!

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Section 2: Plasma activation and chip binding

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Section 3: Coating of the channels

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1. Place the dishes containing the chips in the vacuum jar-bell and turn the pump on. After 1 minute, close the chamber and stop the pump. Wait for 10 to 15 minutes.

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2. Take the dishes under the cell culture hood.

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3. Pipette 10 μ L of fibronectin solution (25 μ g/mL) and add it in each entry holes of the chips as well as on the side of the chip. You should see liquid spreading in the whole structure.

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4. Keep the samples closed outside of the hood in a big plastic Petri dish for 1 hour at room temperature.

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5. After adsorption, under the hood, empty the holes by aspiration. Cast phosphate buffer solution (PBS) in the dish, remove the bubbles in the entry holes and aspirate immediately all the liquid.

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7. After three step of rinsing, empty the whole dish and cast culture medium in it. Half of the samples will have a normal medium while the other half will have a medium containing a special molecule. Put enough medium so that the PDMS is totally covered.

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Section 4: Cell staining, recovery, counting and seeding

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1. Add 5 µL of Hoechst solution in two wells of the plate containing the cells. Replace the plate in the incubator for 30 minutes.

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2. With a 1 mL pipette, aspirate the liquid of one well and do several back and forth movements of the liquid to unstick attached cells. During that step be cautious to not form

bubbles. Put the cell solution in a first tube. Redo this step for the second well and fill a second tube.

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
3. Prepare the Malassez for the counting of cells and deposit 20 μL of a cell solution in one chamber. Calculate the concentration of cells.

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Malassez cell is composed by a grid of 25 rectangles of 20 cases (see scheme 2). In order to count a cell solution, count the amount of cells in 10 rectangles. Multiply this number by 10 and you have the number of cells per μL .

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Scheme2Malassezcell

Abbildung 2 Malassez-Zelle

4. Centrifuge the tubes for 5 minutes at 1200 rpm.

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5. Aspirate the liquid above the cell pellet that has formed at the bottom of the tube. Be careful not to aspirate the pellet itself!

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6. Resuspend the pellet in the correct medium (normal medium for the first pellet and medium containing the Y27632 for the second pellet). Put as much medium as needed to obtain a final concentration of cells of 15 million per mL.

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7. Take samples out of the incubator and lower the level of medium below the level of PDMS and empty the entry holes.

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8. Push 10 mL of cell solution in both entries of one single loading channel, letting the other two entries empty. This step creates a pressure difference which helps to drive the cells in front of the migrating channels.

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9. Check on the microscope that enough cells are in front of the migrating channels and recover the entire chip with medium.

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Now your cells are ready to be imaged on the microscope. This step will be shown and explained directly on site on the day of the experiment. Cells will be observed and recorded overnight.

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Planned timing of steps:

- Chip preparation (45 min)
- Assembly and oven (1h, 20 min for mixing incl.)
- Vacuum and coating (30 min +1h coating: lunch)
- Rinsing, medium, cell staining (20+30+30)
- Cell recovery and counting (30 min)
- Seeding (20 min +30 min for entry)
- Microscopy (1hour)

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Day 2: Data Analysis

Tag 2: Datenanalyse

In this 2nd part of the practical you will develop and discuss representative figures and statistical analysis of your data starting from images taken with the microscope the night before. Therefore, you will follow certain steps along which we have prepared some more questions. Please remember to give the source of your information when you prepare the answers of those questions.

Image analysis will be accomplished by the free software ImageJ, Statistics and Data representation will be performed in Origin. Both programs are available on a PC in the AG Lautenschläger which can be used at 100% at day 2 of the practical and if necessary on additional days in agreement with the group. These programs will be sufficient to analyze manually the measured data. If there is spare time and interest, students are invited to write a small analysis programs to make data analysis more efficient.

We will assist you with the following steps for the first 5 cells. After this, you will repeat the steps yourself until your dataset is large enough to get a statistical answer.

Can we theoretically reach a number of measurements which allows to be sure of our conclusions? How do we deal with this is practically?

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Workflow day 2:

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1. Data transfer from microscope to PC

Datentransfer von Mikroskop zu PC

This step can take up to one hour. In the meantime, discuss the following questions with the person supervising the practical:

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An image is a matrix of values. The intensity measured by each pixel of the camera is expressed as an intensity level compromised in a range depending on the image depth. This depth determines the size of the image once recorded on a computer, and a good balance between file size and intensity resolution has to be found.

How are sizes of objects calibrated in images? What is a kymograph and what is it used for?

Further on, get familiar with ImageJ.

(open, close, save, measure parameters, set scale bar, build movies out of a series of images, generate a kymograph)

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2. Assembly of movies

Filmmontage

Assemble one movie per microscope position (on each position you will have several cells to analyze)

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3. Generation of one kymograph for each movie for following analysis.

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Generate one representative kymograph with a scalebar for your report.

In the kymographs, you will see a different behaviour of movement for each cell. For the sake of simplicity, you will measure the velocity of the cells exhibiting a long straight migration without touching any other cells. Nevertheless, it is important to know how representative of the population those cells are.

Present your results as a table containing for each analyzed kymograph:

- how many cells there are in total
- From this total, how many exhibit a change in direction.
- How many cells did you measure
- From the measured cells, how many have changed direction

What is called persistence in cell migration ?

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4. Statistics

Statistik

Use Origin to measure a mean value with standard error of cell speed for both cell populations. Illustrate the distribution of single cell speed for each population as well as the mean value. Try at least two different plot types and discuss advantages or disadvantages.

Perform a t-test in order to test of differences between both cell populations.

What is a t-test?

Under which conditions a t-test can be performed?

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Interpretation der Daten

Summarize the experiment, your expectations of data and the actual results you obtained. Is the cell speed of the two populations different? If yes or no: what might be a possible explanation? What is the drug/molecule doing within the cell and how can this explain the data you got? Try to find references in Literature to support your results.

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Possible sources to help answering your questions:

Molecular Biology of the Cell by Bruce Alberts *et al.*

www.wikipedia.com (in German AND English)

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Guidelines for report:

Don't forget to answer each question!

Organize the text in different sections.

Don't forget scale bars (for example in kymograph).

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