

## Mouse genotyping with Thermo F170S Phire Tissue Direct PCR Kit (hair follicles)

DNA extraction mix, per sample:

10 $\mu$ L Dilution buffer + 0,5  $\mu$ L DNA Release Additive.

Prepare the required volume of DNA extraction mix. Distribute 10 $\mu$ L in 8-tubes PCR stripes. Place 10-50 hair follicles (wipe the forceps with 70% EtOH paper tissue between sampling) into each tube.

Incubate at room temperature 5 – 30 min. Heat inactivate 96°C 2 min (on PCR cyclor).

Prepare PCR mix (per 10 samples):

50 $\mu$ L Phire Tissue Direct PCR Master Mix (2X)

5 $\mu$ L Forward primer 10 $\mu$ M

5 $\mu$ L Reverse primer 10 $\mu$ M

30 $\mu$ L water (to 90 $\mu$ L)

Distribute 9 $\mu$ L per PCR tube and add 1 $\mu$ L from DNA extraction reaction. Run 33-35 cycles. For annealing temperature calculation use the online tool:

<https://www.thermofisher.com/de/de/home/brands/thermo-scientific/molecular-biology/molecular-biology-learning-center/molecular-biology-resource-library/thermo-scientific-web-tools/tm-calculator.html>

Run/detect/analyse PCR products on 2 – 3% agarose gel.

Alternative protocol (for robust PCR tests, to save 50% of PCR Master Mix (2X)).

PCR mix (per 10 samples):

25 $\mu$ L Phire Tissue Direct PCR Master Mix (2X)

5 $\mu$ L Forward primer 10 $\mu$ M

5 $\mu$ L Reverse primer 10 $\mu$ M


4 $\mu$ L dNTPs (2,5mM each) / 1 $\mu$ L dNTPs (10mM each)

5 $\mu$ L Taq Buffer x 10 (or 10 $\mu$ L Phusion/Q5 buffer x 5)

X $\mu$ L water (to 90 $\mu$ L)



Hair sampling procedure



Example: 200bp (WT) and 178bp (mutant allele, 22bp deletion) fragments of tet3 gene.