

Users at Saarland University can access the MNOVA software. Corresponding campus licences are available on the NMR server for authorised users. The software version 14 can be downloaded for the required operating system from the provider's homepage (<https://mestrelab.com/download/mnova/>). After installation, simply copy the licences received into the "licenses" subfolder in the software installation folder before opening the software for the first time. If required, the plugins for "qNMR" and "Reaction Monitoring" can be installed via the menu item "Advanced Plugins". To activate the licences, the licence server must be contacted. This happens automatically as soon as a connection to the university network is established (VPN connections are also possible). Licences remain valid for 90 days before the licence server has to be contacted again. If you have any questions, please contact our NMR Service Center.